

## HOXA10 antibody

### Product Information


Catalog No.:	FNab10947
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	Homeobox protein Hox-A10
Synonyms:	Homeobox protein Hox-A10 Homeobox protein Hox-1.8 Homeobox protein Hox-1H PL HOXA10 HOX1H
Observed MW:	42-50 kDa
Uniprot ID :	P31260

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis using FNab10947(HOXA10 Antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10947( HOXA10 Antibody) at dilution of 1:1000.