

PTPRF antibody

Product Information

Catalog No.:	FNab10951
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4H3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	protein tyrosine phosphatase, receptor type, F
Synonyms:	Receptor-type tyrosine-protein phosphatase F Leukocyte common antigen related (LAR) PTPRF LAR
Observed MW:	212 kDa
Uniprot ID :	P10586

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10951(PTPRF Antibody) at dilution of 1:1000.