

Phospho-c-Jun(S73) antibody

Product Information


Catalog No.:	FNab10952
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	A synthetic phosphorylated peptide around S73 of human c-JUN
Synonyms:	Transcription factor Jun Activator protein 1 (AP1) Proto-oncogene c-Jun Transcription factor AP-1 subunit Jun V-jun avian sarcoma virus 17 oncogene homolog p39 JUN
Observed MW:	44 kDa
Uniprot ID :	P05412

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain using FNab10952 at dilution of 1:100



PC-14 cells tissue were subjected to SDS PAGE followed by western blot with FNab10952 at dilution of 1:1000