

## PDGFRA antibody

### Product Information


Catalog No.:	FNab10965
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Platelet-derived growth factor receptor alpha
Synonyms:	Platelet-derived growth factor receptor alpha (PDGF-R-alpha, PDGFR-alpha) Alpha platelet-derived growth factor receptor Alpha-type platelet-derived growth factor receptor CD140 antigen-like family member A CD140a antigen Platelet-derived growth factor alpha receptor Platelet-derived growth factor receptor 2 (PDGFR-2) PDGFRA PDGFR2 RHEPDGFRA
Observed MW:	150 kDa
Uniprot ID :	P16234

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB : 1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10965(PDGFRA antibody) at dilution of 1:1000