

THRB antibody

Product Information


Catalog No.:	FNab10968
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Thyroid hormone receptor beta
Synonyms:	Thyroid hormone receptor beta Nuclear receptor subfamily 1 group A member 2 c-erbA-2 c-erbA-beta THRB ERBA2 NR1A2 THR1
Observed MW:	52 kDa
Uniprot ID :	P10828

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human gastric carcinoma tissue using FNab10968 at dilution of 1:100