

LHX2 antibody

Product Information


Catalog No.:	FNab10974
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	LIM homeobox 2
Synonyms:	hLhx2 Homeobox protein LH 2 LH2 LHX2 LIM homeobox 2 LIM homeobox protein 2 LIM/homeobox protein Lhx2
Observed MW:	50-55 kDa
Uniprot ID :	P50458

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



A549 cell lysates were subjected to SDS PAGE followed by western blot with FNab10974(LHX2 antibody) at dilution of 1:500