

RAD51C antibody

Product Information


Catalog No.:	FNab10979
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	DNA repair protein RAD51 homolog 3
Synonyms:	R51H3 RAD51 homolog C RAD51 like protein 2 RAD51C RAD51L2
Observed MW:	42kDa
Uniprot ID :	O43502

Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10979(RAD51C antibody) at dilution of 1:1000