

AQP1 antibody

Product Information


Catalog No.:	FNab10981
Conjugation:	None
Size:	100µg
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4B6
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	Aquaporin-1 Aquaporin-CHIP Urine water channel Water channel protein for red blood cells and kidney proximal tubule AQP1 CHIP28
Uniprot ID :	P29972

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:5000; IHC: 1:50-1:500



human heart tissue were subjected to SDS PAGE followed by western blot with FNab10981(AQP1 antibody) at dilution of 1:2000