

## THY1 antibody

### Product Information


Catalog No.:	FNab10982
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9J3
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Thy-1 membrane glycoprotein
Synonyms:	Thy-1 membrane glycoprotein CDw90 Thy-1 antigen THY1 CD90
Observed MW:	25 kDa
Uniprot ID :	P04216

### Application

Reactivity:	Human, Mouse, Rat, Pig, Rabbit
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:1000 - 1:5000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10982(CD90 antibody) at dilution of 1:2000