

CYP7A1 antibody

Product Information


Catalog No.:	FNab10985
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	cytochrome P450, family 7, subfamily A, polypeptide 1
Synonyms:	Cytochrome P450 7A1 CYPVII CYP7 CP7A Cholesterol 7-alpha-monooxygenase
Observed MW:	58 kDa
Uniprot ID :	P22680

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse liver were subjected to SDS PAGE followed by western blot with FNab10985(CYP7A1 antibody) at dilution of 1:1000