

GLP1R antibody

Product Information


Catalog No.:	FNab10994
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Glucagon-like peptide 1 receptor
Synonyms:	GLP-1-R GLP-1R GLP-1 receptor GLP 1R GLP 1 receptor
Observed MW:	53 kDa
Uniprot ID :	P43220

Application

Reactivity:	Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Rat lung were subjected to SDS PAGE followed by western blot with FNab10994(GLP1R antibody) at dilution of 1:1000