

NLGN1 antibody

Product Information


Catalog No.:	FNab11004
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9J3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Neurologin-1
Synonyms:	KIAA1070 neurologin 1 NLGN1
Observed MW:	96-116 kDa
Uniprot ID :	Q8N2Q7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:5000
Image:	



Mouse brain tissue were subjected to SDS
PAGE followed by western blot with
FNab11004 at dilution of 1:2000