

PLIN1 antibody

Product Information


Catalog No.:	FNab11009
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Perilipin 1
Synonyms:	PLIN PLIN1 PERI Perilipin 1 Lipid droplet-associated protein
Observed MW:	60-70 kDa
Uniprot ID :	O60240

Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab11009 at dilution of 1:1000