

GPER1 antibody

Product Information


Catalog No.:	FNab11011
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	G-protein coupled estrogen receptor 1
Synonyms:	mER CEPR GPER DRY12 FEG-1 GPR30 LERGU LyGPR CMKRL2 LERGU2 GPCR-Br GPER1
Observed MW:	42-52 kDa
Uniprot ID :	Q99527

Application

Reactivity:	Human, Mouse
Tested Application:	WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:50 - 1:200
Image:	



Immunofluorescent analysis of NIH/3T3 cells using FNab11011 at dilution of 1:50



Mouse heart were subjected to SDS PAGE followed by western blot with FNab11011 at dilution of 1:1000