

Phospho-DRP1(S637) antibody

Product Information


Catalog No.:	FNab11012
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	peptide
Synonyms:	Dynamin-1-like protein Dnm1p/Vps1p-like protein (DVLP) Dynamin family member proline-rich carboxyl-terminal domain less (Dymple) Dynamin-like protein Dynamin-like protein 4 Dynamin-like protein IV (HdynIV) Dynamin-related protein 1 DNM1L DLP1 DRP1
Observed MW:	82 kDa
Uniprot ID :	O00429

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab11012 at dilution of 1:1000