

## Phospho-p53 (S392) antibody

### Product Information


Catalog No.:	FNab11015
Conjugation:	None
Size:	100µg
Form:	liquid
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	Cellular tumor antigen p53 Antigen NY-CO-13 Phosphoprotein p53 Tumor suppressor p53 TP53 P53
Uniprot ID :	P04637

### Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



+   -   phosphatase inhibitor  
-   +    $\lambda$  phosphatase

phosphatase inhibitor treated and  $\lambda$   
phosphatase treated HEK293 cells were  
subjected to SDS PAGE followed by western  
blot with FNab11015 at dilution of 1:1000