
SLC7A5 antibody

Product Information


Catalog No.:	FNab11016
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6F8
IsoType:	IgG2b
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A synthetic phosphorylated peptide around S392 of human p53.
Synonyms:	Cellular tumor antigen p53 Antigen NY-CO-13 Phosphoprotein p53 Tumor suppressor p53 TP53 P53
Observed MW:	40-50 kDa
Uniprot ID :	Q01650

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab11016 at dilution of 1:1000