

## SFRP1 antibody

### Product Information


Catalog No.:	FNab11019
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Secreted frizzled-related protein 1
Synonyms:	FRP FRP 1 FRP1 FrzA SARP 2 SARP2 sFRP 1 SFRP1
Observed MW:	30-35 kDa
Uniprot ID :	Q8N474

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab11019 at dilution of 1:1000