

KRT5 antibody

Product Information

Catalog No.:	FNab11020
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells, which are classified into two major sequence types. This antibody is a pan-keratin antibody.

Immunogen information


Immunogen:	keratin 5
Synonyms:	Keratin, type II cytoskeletal 5 58 kDa cytokeratin Cytokeratin-5 (CK-5) Keratin-5 (K5) Type-II keratin Kb5 KRT5
Observed MW:	60 kDa
Uniprot ID :	P13647

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab11020 at dilution of 1:100



A375 cells were subjected to SDS PAGE followed by western blot with FNab11020 at dilution of 1:1000