

CD133 antibody

Product Information


Catalog No.:	FNab11021
Conjugation:	None
Size:	100µg
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2H6
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Prominin-1 Antigen AC133 Prominin-like protein 1 PROM1 PROML1 CD133
Uniprot ID :	O43490
Observed MW:	115 kDa-120 kDa

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab11021 at dilution of 1:200



HT29 cells were subjected to SDS PAGE followed by western blot with FNab11021 at dilution of 1:1000