

Phospho-LRP6(S1490) antibody

Product Information


Catalog No.:	FNab11029
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Low-density lipoprotein receptor-related protein 6
Synonyms:	Low-density lipoprotein receptor-related protein 6 LRP6 LRP-6
Observed MW:	180 kDa
Uniprot ID :	O75581

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:	WB: 1:500 - 1:1000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab11029 at dilution of 1:500. Hela cells were treated by Calyculin A (100 nM) at 37°C for 30 minutes after serum-starvation overnight.