

Phospho-RPTOR(S792) antibody

Product Information


Catalog No.:	FNab11030
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	raptor
Synonyms:	Regulatory-associated protein of mTOR (Raptor) p150 target of rapamycin (TOR)-scaffold protein RPTOR KIAA1303 RAPTOR
Observed MW:	150 kDa
Uniprot ID :	Q8N122

Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



293T cells were subjected to SDS PAGE followed by western blot with FNab11030 at dilution of 1:1000; 293T cells were treated by Oligomycin (0. 5 uM) at 37°C for 30min after serum-starvation overnight.