

---

## IL-15 antibody

### Product Information


Catalog No.:	FNab11032
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	full length protein
Synonyms:	IL15 IL-15
Observed MW:	15 kDa
Uniprot ID :	P40933

### Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Mouse heart were subjected to SDS PAGE followed by western blot with FNab11032 at dilution of 1:1000