

PRKN antibody

Product Information


Catalog No.:	FNab11037
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7J6
IsoType:	IgG2b
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	parkin
Synonyms:	E3 ubiquitin-protein ligase parkin (Parkin) Parkin RBR E3 ubiquitin-protein ligase Parkinson juvenile disease protein 2 (Parkinson disease protein 2) PRKN PARK2
Observed MW:	52 kDa
Uniprot ID :	O60260

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse brain using FNab11037 Antibody at dilution of 1:100



U251 cells were subjected to SDS PAGE followed by western blot with FNab11037 at dilution of 1:1000