

LDHAL6B antibody

Product Information


Catalog No.:	FNab11041
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	L-lactate dehydrogenase A-like 6B
Synonyms:	LDHL DH6B LDHAL6 LDHAL6B
Observed MW:	42 kDa
Uniprot ID :	Q9BYZ2

Application

Reactivity:	Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:50-1:200
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab11041 at dilution of 1:1000