

Substance P antibody

Product Information


| | |
|---------------|--|
| Catalog No.: | FNab12001 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Immunogen information

| | |
|--------------|---|
| Immunogen: | KLH conjugated synthetic peptide derived from human Substance P |
| Synonyms: | Protachykinin-1 PPT Substance P Neurokinin A (NKA) Alternative names: Neuromedin L Substance K Neuropeptide K (NPK) Neuropeptide gamma C-terminal-flanking peptide TAC1 NKA NKNA TAC2 |
| Observed MW: | None |
| Uniprot ID : | P20366 |

Application

| | |
|-----------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, IHC |
| Recommended dilution: | IHC: 1:50 - 1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded rat brain tissue using FNab12001 at dilution of 1:100