

---

## Purified Anti-Mouse CD3E (145-2C11)

### Product Information

Catalog No.:	FNab30002
Size:	100µg
Form:	liquid
Purity:	≥95% as determined by SDS-PAGE
Host:	Armenian Hamster
Clonality:	monoclonal
Clone ID:	145-2C11
IsoType:	IgG
Storage:	PBS, pH 7.2. Store at 2-8° C valid for 12 months, avoid freeze/thaw cycles.

### Immunogen information

Synonyms:	T-cell surface glycoprotein CD3 epsilon chain  CD3E  T-cell surface antigen T3/Leu-4 epsilon chain  CD3e  CD3E  T3E
Uniprot ID :	P22646

### Application

Reactivity:	Mouse
Tested Application:	FC
Recommended dilution:	2 µg/mL(1×10 <sup>5</sup> -5×10 <sup>5</sup> cells)