
Purified Anti-Mouse CD28(37.51)

Product Information

Catalog No.:	FNab30019
Size:	100µg
Form:	liquid
Purity:	≥95% as determined by SDS-PAGE
Host:	Syrian Hamster
Clonality:	monoclonal
Clone ID:	37.51
IsoType:	IgG
Storage:	PBS, pH 7.2. Store at 2-8° C valid for 12 months, avoid freeze/thaw cycles.

Immunogen information

Synonyms:	T-cell-specific surface glycoprotein CD28 Cd28 CD28
Uniprot ID :	P31041

Application

Reactivity:	Mouse
Tested Application:	FC
Recommended dilution:	2 µg/mL(1×10^5 - 5×10^5 cells)