
Purified Anti-Mouse IFN γ (XMG1.2)

Product Information

Catalog No.:	FNab30074
Conjugation:	None
Size:	100 μ g
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	XMG1.2
IsoType:	IgG1, kappa
Storage:	PBS, pH 7.2. Store at 2-8°C valid for 12 months, avoid freeze/thaw cycles.

Immunogen information

Synonyms:	Interferon gamma IFNG IFN-gamma IFN γ IFN- γ
Uniprot ID :	P01580

Application

Reactivity:	Mouse
Tested Application:	FC
Recommended dilution:	2 μ g/mL(1×10^5 - 5×10^5 cells)