

FineTest®Violet 500 Anti-Mouse MHCII Antibody(M5/114)

Product Information

Catalog No.:	FV500-30124
Conjugation:	FineTest®Violet 500
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	M5/114
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	H2-Ab1 Eb1 I-E beta MHC class II MHC class II MHC class II H2-IA-beta-psi Major histocompatibility protein class II beta chain
Uniprot ID :	P14483 O78196

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.