

FineTest®Violet 540 Anti-Mouse CD8A (53-6.7)

Product Information

Catalog No.:	FV540-30003
Conjugation:	FineTest®Violet 540
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	53-6.7
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	T-cell surface glycoprotein CD8 alpha chain CD8A T-lymphocyte differentiation antigen T8/Leu-2 MAL
Uniprot ID :	P01731

Application

Reactivity:	Mouse
Tested Application:	IF, FC
Recommended dilution:	IF/ICC : 1:250-1:1000; FC:0.1-1 µg/10 ⁶ cells in 100 µL volume.
Image:	