

HRP-GAPDH monoclonal antibody

Product Information

Catalog No.:	HRP-03343
Size:	100µg
Form:	liquid
Purification:	Antigen affinity purification
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4A2
IsoType:	IgG1
Storage:	PBS with 0.05% Proclin300, 0.5% BSA and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	glyceraldehyde-3-phosphate dehydrogenase
Synonyms:	G3PD, GAPD, GAPDH
Observed MW:	36 kDa
Uniprot ID :	P04406

Application

Reactivity:	Human, Mouse, Rat, Zebrafish
Tested Application:	WB
Recommended dilution: WB:	1:5000-1:50000