

PerCP-Cy5.5 Anti-Human CD8A (OKT-8)

Product Information

Catalog No.:	PC55-30006
Conjugation:	PerCP-Cy5.5
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	OKT-8
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	T-cell surface glycoprotein CD8 alpha chain CD8A T-lymphocyte differentiation antigen T8/Leu-2 MAL
Uniprot ID :	P01732

Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Human peripheral blood lymphocytes are stained with PerCP/Cy5.5 Anti-Human CD8a Antibody (red histogram). Unstained lymphocytes (empty black histogram) are used as control.