

PE-Cy7 Anti-Mouse CD80(16-10A1)

Product Information

Catalog No.:	PC7-30024
Conjugation:	PE-Cy7
Size:	50T/100T
Form:	liquid
Host:	Armenian Hamster
Clonality:	monoclonal
Clone ID:	16-10A1
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	T-lymphocyte activation antigen CD80 Cd80 Activation B7-1 antigen B7
Uniprot ID :	Q00609


Application

Reactivity: Mouse

Tested Application: FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



RAW264.7 cells are stained with PE/Cy7 Anti-Mouse CD80 Antibody (filled red histogram) or PE/Cy7 Armenian Hamster IgG Isotype Control (empty black histogram).

