

## PE-Cy7 Anti-Human CD117(104D2)

### Product Information

Catalog No.:	PC7-30047
Conjugation:	PE-Cy7
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	104D2
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Synonyms:	Mast/stem cell growth factor receptor Kit Kit SCFR Proto-oncogene c-Kit Tyrosine-protein kinase Kit CD117
Uniprot ID :	P10721

### Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



HEL cells are stained with PE/Cy7 Anti-Human CD117 Antibody (filled red histogram). Unstained cells (empty black histogram) are used as control.