

---

## PerCP/Cy5.5 Anti-Mouse CD49b(pan-NK cells) Antibody(DX5)

### Product Information

Catalog No.:	PCP55-30097
Conjugation:	PerCP/Cy5.5
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	DX5
IsoType:	Rat IgM, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	CD49 antigen-like family member B CD49b Collagen receptor GPIa Integrin alpha-2 Platelet membrane glycoprotein Ia VLA- 2 subunit alpha
Uniprot ID :	Q62469

### Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.