

PerCP-Cy5.5 Armenian Hamster IgG Isotype Control(PIP)

Product Information

Catalog No.:	PCP55-30104
Conjugation:	PerCP-Cy5.5
Size:	20T/100T
Form:	liquid
Host:	Armenian Hamster
Clonality:	None
Clone ID:	PIP
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	None
Uniprot ID :	None

Application

Reactivity:	None
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: