

---

## PE Anti-Mouse CD4(GK1.5)

### Product Information

Catalog No.:	PE-30001
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	GK1.5
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	T-cell surface glycoprotein CD4, T-cell surface antigen T4, Leu-3, L3T4, T4
Uniprot ID :	P06332

### Application

Reactivity:	Mouse
Tested Application:	FC
Recommended dilution: Volume per test:	5 $\mu$ L