
PE Anti-Mouse CD3E(145-2C11)

Product Information

| | |
|--------------|----------------------------------------------------------------------------------------------------------|
| Catalog No.: | PE-30002 |
| Conjugation: | PE |
| Size: | 50T/100T |
| Form: | liquid |
| Host: | Armenian Hamster |
| Clonality: | monoclonal |
| Clone ID: | 145-2C11 |
| IsoType: | IgG |
| Storage: | PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Immunogen information

| | |
|--------------|---------------------------------------------------------------------------------------------------------------------|
| Synonyms: | T-cell surface glycoprotein CD3 epsilon chain, CD3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T3E |
| Uniprot ID : | P22646 |

Application

| | |
|----------------------------------------|-------|
| Reactivity: | Mouse |
| Tested Application: | FC |
| Recommended dilution: Volume per test: | 5 µ L |