

---

## PE Anti-Human CD4 (RPA-T4)

### Product Information

Catalog No.:	PE-30005
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	RPA-T4
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	T-cell surface glycoprotein CD4, T-cell surface antigen T4, Leu-3, L3T4, T4
Uniprot ID :	P01730

### Application

Reactivity:	Human
Tested Application:	FC
Recommended dilution: Volume per test:	5 µ L