

---

## PE Anti-Human CD8A (OKT-8)

### Product Information

Catalog No.:	PE-30006
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	OKT-8
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	T-cell surface glycoprotein CD8 alpha chain, CD8A, T-lymphocyte differentiation antigen T8/Leu-2, MAL
Uniprot ID :	P01732

### Application

Reactivity:	Human
Tested Application:	FC
Recommended dilution:	Volume per test: 5 µ L