

PE Anti-Mouse CD86(GL-1)

Product Information

Catalog No.:	PE-30025
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	GL-1
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	T-lymphocyte activation antigen CD86 Cd86 Activation B7-2 antigen Early T-cell costimulatory molecule 1 ETC-1
Uniprot ID :	P42082

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



LPS-stimulated (3 days) KM mouse splenocytes are stained with PE Anti-Mouse CD86 Antibody (red histogram). Unstained splenocytes (empty black histogram) are used as control.