

PE Anti-Mouse NK1.1(PK136)

Product Information

Catalog No.:	PE-30029
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	PK136
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Killer cell lectin-like receptor subfamily B member 1C Klrb1c CD161 antigen-like family member C Ly-55c CD161/NK1.1 NKR-P1.9 NKR- P1C NKR-P1 40 CD161c
Uniprot ID :	P27814/P27812/Q99JB4

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



C57BL/6 murine splenocytes are stained with PE Anti-Mouse NK1.1 Antibody and APC Anti-Mouse CD49b Antibody.