

PE Anti-Human/Mouse/Rat CD47(MIAP410)

Product Information

Catalog No.:	PE-30031
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	MIAP410
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	Leukocyte surface antigen CD47 Cd47 Integrin-associated protein IAP
Uniprot ID :	Q08722/Q61735/P97829


Application

Reactivity: Human,Mouse,Rat

Tested Application: FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



C57BL/6 murine splenocytes are stained with PE Anti-Human/Mouse/Rat CD47 Antibody (filled red histogram). Unstained splenocytes (empty black histogram) are used as control.

