

PE Anti-Human CD117(104D2)

Product Information

| | |
|--------------|--|
| Catalog No.: | PE-30047 |
| Conjugation: | PE |
| Size: | 50T/100T |
| Form: | liquid |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 104D2 |
| IsoType: | IgG1, kappa |
| Storage: | PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Immunogen information


| | |
|--------------|---|
| Synonyms: | Mast/stem cell growth factor receptor Kit Kit SCFR Proto-oncogene c-Kit Tyrosine-protein kinase Kit CD117 |
| Uniprot ID : | P10721 |

Application

| | |
|---------------------|-------|
| Reactivity: | Human |
| Tested Application: | FC |

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



HEL cells are stained with PE Anti-Human CD117 Antibody (filled red histogram). Unstained HEL cells (empty black histogram) are used as control.