

## PE Anti-Human CD122(TU27)

### Product Information

Catalog No.:	PE-30051
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	TU27
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Synonyms:	Interleukin-2 receptor subunit beta IL-2RB p75 CD122 IL-2R $\beta$
Uniprot ID :	P14784

### Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5 $\mu$ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5  $\mu$ L of antibody per test (million cells in 100  $\mu$ L staining volume or per 100  $\mu$ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Human peripheral blood lymphocytes are stained with APC Anti-Human CD56 Antibody and PE Anti-Human CD122 Antibody(Left). Lymphocytes are stained with APC Anti-Human CD56 Antibody and PE Mouse IgG1,  $\kappa$  Isotype Control (Right).

