

PE Anti-Mouse IFN γ (XMG1.2)

Product Information

Catalog No.:	PE-30074
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	XMG1.2
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Interferon gamma IFNG IFN-gamma IFN γ IFN- γ
Uniprot ID :	P01580

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5 μ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



HEK293T cells transiently transfected with pcDNA3.1 plasmid encoding Mouse IFN- γ gene are stained with PE Anti-Mouse IFN- γ Antibody (red histogram) or PE Rat IgG1, κ Isotype Control (black histogram).