

PE Anti-Human/Mouse CD44 Antibody(IM7)

Product Information

Catalog No.:	PE-30083
Conjugation:	PE
Size:	20T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	IM7
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	CD44 antigen CD44 CDw44 Epican Phagocytic glycoprotein 1 PGP-1 Phagocytic glycoprotein I PGP-I CD44 LHR MDU2 MDU3 MIC4
Uniprot ID :	P15379/P16070

Application

Reactivity:	Human/Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



KM murine splenocytes are stained with PE Anti-Human/Mouse CD44 Antibody (red histogram). Unstained splenocytes (black histogram) are used as control.