

## PE Anti-Mouse CD11C Antibody(N418)

### Product Information

Catalog No.:	PE-30085
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Armenian Hamster
Clonality:	monoclonal
Clone ID:	N418
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	Integrin alpha-X Itgax CD11 antigen-like family member C Leukocyte adhesion receptor p150+95 CD11c
Uniprot ID :	Q9QXH4

### Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: