

PE Anti-Human CD197 Antibody(G043H7)

Product Information

Catalog No.:	PE-30090
Conjugation:	PE
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	G043H7
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	C-C chemokine receptor type 7 (C-C CKR-7, CC-CKR-7, CCR-7) BLR2 CDw197 Epstein-Barr virus-induced G-protein coupled receptor 1 (EBI1, EBV-induced G-protein coupled receptor 1) MIP-3 beta receptor CCR7 CMKBR7 EBI1 EVI1
Uniprot ID :	P32248

Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Human peripheral blood lymphocytes are stained with APC Anti-Human CD3 Antibody and PE Anti-Human CD197 Antibody (Left). Lymphocytes are stained with APC Anti-Human CD3 Antibody and PE Mouse IgG2a, κ Isotype Control (Right).